

signia

Life sounds brilliant.

Hear the colors of sound.

The new Prompt hearing aids.
For the important things in life.



Hearing
Systems

SIEMENS

signia.to/prompt

The new Prompt.

Whatever your wearing preference, the new Prompt hearing aids offer a range of solutions for almost all types of hearing loss. Equipped with advanced technology such as state-of-the-art feedback cancellation, Prompt provides great sound. Three behind-the-ear or BTE models offer outstanding wearing comfort together with ease of use. With their ingenious, interchangeable, super-soft Click sleeves that simply click on the hearing aid, the in-the-ear or ITE models fit beautifully and in seconds. With their range of selectable programs, Prompt hearing aids help you in the most diverse hearing situations.



Some things are made for sharing.

From a whispered compliment to the sound of delighted children laughing or summer birdsongs in the garden, life is full of moments meant for sharing. With Prompt hearing aids there to help, you can enjoy every moment with complete confidence and be a part of it.



Great for every conversation, whoever it's with.



At home, whether it's quiet or the kids are running around.



In crowded, very noisy environments, such as restaurants or marketplaces.

Prompt BTEs



Prompt S Prompt P Prompt SP

Prompt Click ITEs



Prompt Click CIC Prompt Click ITC

Two different models are available:



A snug fit:

Prompt behind-the-ear models sit comfortably behind the ear.



Absolutely discreet:

Prompt in-the-ear models are worn directly and discreetly in the ear.

Discover Prompt hearing aids in detail:

Prompt BTEs.

A versatile range for optimum amplification.

Prompt Click ITEs.

One click is all it takes for a perfect fit.

UNDERSTAND CONVERSATION

- 1 Focus on what is really important.**
Prompt BTEs' directional microphones pick up the sounds that will be amplified. They focus on sound coming from in front of you while reducing noises from other directions.

EASY HANDLING

- 2 All options at the tip of your finger.**
The push button and/or the rocker switch on Prompt BTEs let you easily and quickly change the program if your listening situation changes or adjust the volume.




- 3 Rugged, robust and reliable.**
Prompt BTEs' housings are trialed and tested in the most challenging environments and are resistant to moisture, sweat, dust and dirt to keep out the rigors of everyday life.

SOUND TRANSMISSION

- 4 Your powerful connection to life.**
Prompt P and SP models come with an earhook for high amplification. The lightweight earmold option provides a secure fit that directs sound straight into your ears.

QUICK-CHANGE ARTIST

- 5 Back on track in 3 simple steps.**
The battery door on Prompt BTEs can be effortlessly opened and closed so you can change the standard battery quickly and easily as the following steps show.



HEAR THE SPEAKER CLEARLY

- 6 Listen closely, even at a distance.**
In places equipped with an induction loop, the Telecoil or T-Coil function in Prompt P and SP models helps you understand speakers more clearly. Simply select the T-program wherever you see the following sign:



- 7 Personalizing whom you hear.**
With Prompt P and SP, you can change the battery housing for an optional audio shoe. When there's a lot of ambient sound, the speaker's voice is streamed directly to your ears via a microphone.



Colors for Prompt BTEs:

All Prompt BTE models are available in three discreet colors so they are almost unnoticeable behind your ear.



HEAR BETTER

- 1 Louder and clearer sound.**
Sounds are transmitted via the microphone and amplified for clearer hearing.

A MORE NATURAL FEELING


- 2 A breath of fresh air.**
The vent in the Click sleeve provides air circulation to the ear to give you the feeling of a normal open ear canal.

Colors for Prompt Click ITEs:

Available in natural skin colors of the ear, Prompt Click ITEs are super discreet and almost invisible in the ear.



Prompt P



Prompt Click CIC

SIMPLE TO INSERT AND REMOVE

- 3 Quick, convenient fitting.**
Easily interchangeable, Click sleeves are made from super-soft silicone which quickly and simply clicks on to the device for a superbly comfortable and swift fit. Available in three different sizes, there's a Click sleeve for almost every ear.
- 4 Handy removal cord.**
Discreet yet convenient, the removal cord lets you easily remove the hearing aids from your ear.

CHANGING BATTERIES

- 5 Battery housing.**
Easy to open so you can quickly and simply exchange the standard battery.



Your Hearing Care Professional will help you decide which hearing aid is the right one for you, depending on your hearing loss, personal requirements and wishes.

The information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases and are subject to change without prior notice.

Legal Manufacturer

Signia GmbH
Henri-Dunant-Strasse 100
91058 Erlangen
Germany

Signia GmbH is a Trademark Licensee of Siemens AG.

Order No. 02928-99C1-7600 | Printed in Germany
© 07.2016 Signia GmbH

Find us on:

signia-hearing.com

facebook.com/signiahearing

youtube.com/signiahearing